



PRODUCT CATALOG

Intelligent Innovation

ONCE MEDICAL CO.,LTD

923/70P. SOI EAKACHAI30 EAKACHAI Rd. MAHACHAI
MUANG SAMUTSAKHON THAILAND 74000
TEL:+66(0) 34 422 449 FAX:+66(0) 34 413514 (AUTO)
WWW.ONCEMEDICALDEVICE.COM



Contents

Introductions *page 2*

ONCE Syringe Technical Information *page 3*

ONCE Disposable Syringe *page 4*

ONCE Safety Syringe *page 5*

Needle-Stick injuries fact *page 6*

ONCE Painless *page 7*

ONCE Concept syringe *page 8*

ONCE Hypodermic needle *page 9*



Introduction

We believe a great company culture is the key to our success. As we've grown over the years, we felt it was important to try and distill our core values and what we feel are important to us as a business, the way we work and the way we think;

ONCE was established in Thailand since year 1999, with the local government support in all aspect of certification and distribution. It has advanced to be one of the largest medical syringes , Hypodermic needles, Safety syringes include blood collection tubes supplier in Thailand and oversea since year 2010. ONCE has a network of professional distributors who are highly regarded in the market place and who provide added value to the ONCE product and services.

ONCE Brand Positioning

ONCE are synonym as ONE, New technology, Creative, Excellence



Our Core Brand Values

Trustworthy

We strive our products are of the highest quality and our brand trustable by our costumer with state-of-the-art technology, fully automated processes and continuous research and development that guarantee the quality of our products are at the highest range.

Innovation

In our effort to offer more, we embrace the innovation of our processes and products to meet the ever changing needs of customer. It is our goal to ensure satisfaction with new and improved products at affordable prices.

Management Team

The quality of a company's management team is one of the best predictors of success. We always believe words doesn't symbolized everything, but with our proven track records in winning the main hospital tender and coverage network within the past year should be enough showing our commitment in this industry.

ONCE™ Syringe technical information

All ONCE™ Disposable syringe consist all the aspect of international materials guideline , the quality which been certified by ISO13485, CE0197 They have features of tips like as central luer slip/ central luer lock/ off-central luer slip, pyrogen free, sterilized by EO gas, Disposable use, Valid for 5 years.



ONCE Luer-Lok™ Tip

Generally used for injections requiring a secure connection of the syringe to another device. The tip is threaded for a 'locking' fit, and is compatible with a variety of needles, catheters, and other devices.

ONCE Luer-Lok™ Tip Sizes:
1 ML
3 ML
5 ML
10 ML
20 ML
30 ML
60 ML



Luer Slip Tip

A friction-fit connection that requires the clinician to insert the tip of the syringe into the needle hub or other attaching device in a push-and-twist manner. This will ensure a connection that is less likely to detach. Simply sliding the attaching device onto the syringe tip will not ensure a secure fitting.


Slip Tip Sizes:
1 ML
3 ML
5 ML
10 ML
20 ML
30 ML
60 ML



Eccentric Luer Slip Tip

Allows for work requiring closer proximity to the skin. Generally used for venipunctures and aspiration of fluids. (Also see luer slip instructions above).


Eccentric Tip Sizes:
10 ML
20 ML
60 ML



Catheter Tip

Used for flushing (cleaning) catheters, gastrostomy tubes and other devices. Insert catheter tip securely into catheter or gastrostomy tube. If leakage occurs, refer to your facility's guidelines.

Catheter Tip Sizes:
2 oz. only



Above shown the common label of ONCE disposable syringe which often use by healthcare staff.

ONCE™ Disposable syringe

Type of tube maternal: polypropylene (Clear type), incorporate the Latex-Free tip of plunger which have stop-line for when pull the plunger at the end of bottom in easy use; emphasize indicator and have the ONCE™ features as:



ONCE tip of plunger design assist to reduce dead space of tip of syringe in order to precision dosage



- Removable Plunger
- Protect reuse of syringe
 - Decrease medical waste product



- Decrease Inner diameter tip of syringe improve precision dosage
- Later free rubber with low resistance edge which is easy to push and pull the plunger

Syringe 0.5ml, 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml

PACKING					ITEM NO.
PRODUCT	QUANTITIES				
	PCS/INNER BOX	PCS/CARTON	REF.CARTON SIZE	WEIGHT/CARTON	
0.5 ML	100 PCS	3200 PCS	572x440x490 mm	19000 grams	ONCSS-005
1 ML	100 PCS	2400 PCS	572x440x490 mm	10700 grams	ONCSS-01
2 ML	100 PCS	2400 PCS	430x425x488 mm	14500 grams	ONCSS-02
3 ML	100 PCS	2400 PCS	430x425x488 mm	15700 grams	ONCSS-03
5 ML	100 PCS	1200 PCS	580x332x420 mm	13450 grams	ONCSS-05
10 ML	100 PCS	1200 PCS	574x425x490 mm	14500 grams	ONCSS-10
20 ML	50 PCS	800 PCS	561x415x535 mm	13750 grams	ONCSS-20
30 ML	50 PCS	480 PCS	564x412x522 mm	10400 grams	ONCSS-30
50 ML	25 PCS	400 PCS	564x412x522 mm	10560 grams	ONCSS-50
60 ML	25 PCS	400 PCS	564x412x522 mm	10560 grams	ONCSS-60

ONCE™ Safety syringe

Technology innovation of safety syringe for protect and reduce incidence of NEEDLESTICK INJURY



3 ml ONCE™ Safety Syringe

Type of Safety syringe

- ONCE™ SAFETY SYRINGE 1 ML WITH NEEDLE NO. 27G x 1/2" CODE ONSS-0127/2
- ONCE™ SAFETY SYRINGE 3 ML WITH NEEDLE NO. 23G x 1 1/4" CODE ONSS-03231/4
- ONCE™ SAFETY SYRINGE 5 ML WITH NEEDLE NO. 22G x 1 1/4" CODE ONSS-05221/4




1. Commonly use as conventional syringe with needle.



2. After use, push the plunger into the top of syringe.



3. Pull the plunger down which needle-keeping mechanism worked automatically.



4. Broke the plunger and needle do safety-lock in syringe tube.

ONCE™ Safety Blunt needle



- ONCE™ Safety Disposable Blunt needle 18G x 1 1/4" CODE ONBN-18G

Needle-stick injuries fact, further information please visit our site: WWW.ONCEMEDICALDEVICE.COM

Needlestick injuries fact







- An estimated 600,000 to 800,000 needlestick and other percutaneous injuries among U.S. health care workers occur annually.¹
- At an average U.S. hospital, workers incur approximately 30 needlestick injuries per 100 beds per year.¹
- Studies show that nurses sustain the majority of these injuries and that as many as onethird of all sharps injuries have been reported to be related to the disposal process.²
- One out of every seven U.S. healthcare workers is accidentally stuck by a contaminated sharp every year.³
- It is believed that only one out of three needlesticks are even reported.³
- As a result of sharps injuries, more than three-dozen U.S. healthcare workers a year contract HIV, two thousand workers a year become infected with hepatitis C, and 400 contract hepatitis B.³
- There are more than 100,000 needlestick injuries in UK hospitals each year.³
- Some 10 million needles are used every year in the NHS* in Scotland with around 4,000 needlestick injuries, two -thirds of which are suffered by nurses.³
- Thousands of Australian hospital workers suffer needle-stick injuries each year and the incidence is increasing. Doctors, nurses and scientists warn that the problem is far greater than statistics reveal, claiming that many injuries go unreported.³
- Victorian Work Cover Authority figures show that the number of claims for needle-stick injuries in Australia jumped 12.6 per cent between 1996-97 and 1997-98.³
- **Needlestick injuries are preventable.** Over 80% of needlestick injuries could be prevented with the use of safer needle devices.⁴
- The US laws are likely to be replicated in Japan, Singapore and possibly other countries including the UK.⁵
- 20 or more blood borne pathogens, including **Human Immunodeficiency Virus (HIV)**, **Hepatitis B**, **Hepatitis C**, can be transmitted via contaminated needlestick injuries.

อ้างอิงจาก* National Audit Office study Health & Safety in NHS Acute Hospital Trusts in England

1. Prevention Needlestick Injuries in Health Care Settings, NIOSH Alert, 11/1999.
2. Occupational Safety and Health Administration (OSHA).
3. International Sharps Injury Prevention Society (ISIPS).
4. Nursing Facts: Needlestick Injury, the American Nurses Association 2002.
5. Pharmacy Today; 1 March 2002; Retractable Needle Syringe.

ONCE™ Painless

Because of Painful is main cause to refuse Needle medication (injection) like as Vaccination and Anti-Biotic medication which do so important to protect and relief the symptom. Tiny needle and Triple Bevel cut sharpness help to decrease puncture and reduce pain when penetrate into the skin. ONCE™ Painless designed Wall of Needle to Thin Stainless steel which assist healthcare staff to push/ pull the plunger easily; shorten needle is proper to penetrate into subcutaneous skin vertically when release medication via needle into the muscle it could reduce pain and swelling.

	Product code	Feature	Syringe size	Needle gauge	Needle length	Scale per unit
	ONIS-0503008	Insulin	0.5 ml	30	8 mm.	1:1
	ONIS-0503108	Insulin	0.5 ml	31	8 mm.	1:1
	ONIS-0503208	Insulin	0.5 ml	32	8 mm.	1:1
	ONIS-1003008	Insulin	1 ml	30	8 mm.	1:2
	ONIS-1003108	Insulin	1 ml	31	8 mm.	1:2
	ONIS-1003208	Insulin	1 ml	32	8 mm.	1:2

ONCE Painless design Tiny needle with thin stainless wall.







Triple Bevel Cut sharpness.

Frontal Needle Cap could place up to easily keep the needle.



Protective Plunger - Cap assist to protect contamination when packing broken.















	Product code	Feature	Syringe size	Needle gauge	Needle length	Scale per ml
	ONRS-0125/2	Medication (M)	1 ml	25	½ inch	1:0.01
	ONRS-0127/2W	M/ Vaccination	1 ml	27	½ inch	1:0.01
	ONRS-0127/2B	M/ Vaccination	1 ml	27	½ inch	1:0.01
	ONRS-0129/2	Vaccination	1 ml	29	½ inch	1:0.01









- Hypodermic needle 30G; 8, 13 (½) mm. CODE ONHN-30G/2
- Hypodermic needle 31G; 8, 13 (½) mm. CODE ONHN-31G/2
- Hypodermic needle 32G; 8, 13 (½) mm. CODE ONHN-32G/2

ONCE™ Concept syringe







For Intravenous IV, Intramuscular IM, Emergency medicine

	Product code	Syringe size	Needle gauge	Needle length	Scale per ml
	ONRS-0123	1 ml	23	1 inch	1:0.01
	ONRS-01231/4	1 ml	23	1¼ inch	1:0.01
	ONRS-0124	1 ml	24	1 inch	1:0.01
	ONRS-01241/4	1 ml	24	1¼ inch	1:0.01
	ONRS-0125	1 ml	25	1 inch	1:0.01
	ONRS-01251/4	1 ml	25	1¼ inch	1:0.01
	ONRS-0127W	1 ml	27	1 inch	1:0.01
	ONRS-0127B	1 ml	27	1 inch	1:0.01
	ONRS-0323	3 ml	23	1 inch	1:0.1
	ONRS-03231/4	3 ml	23	1¼ inch	1:0.1
	ONRS-0324	3 ml	24	1 inch	1:0.1
	ONRS-03241/4	3 ml	24	1¼ inch	1:0.1
	ONRS-0325	3 ml	25	1 inch	1:0.1
	ONRS-03251/4	3 ml	25	1¼ inch	1:0.1

For Draw the medicine

	Product code	Syringe size	Needle gauge	Needle length	Scale per ml
	ONRS-0118	1 ml	18	1 inch	1:0.01
	ONRS-01181/4	1 ml	18	1¼ inch	1:0.01
	ONRS-0318	3 ml	18	1 inch	1:0.1
	ONRS-03181/4	3 ml	18	1¼ inch	1:0.1
	ONRS-0518	5 ml	18	1 inch	1:0.2
	ONRS-05181/4	5 ml	18	1¼ inch	1:0.2

For Blood sampling

	Product code	Syringe size	Needle gauge	Needle length	Scale per ml
	ONRS-0521	5 ml	21	1 inch	1:0.2
	ONRS-05211/2	5 ml	21	1½ inch	1:0.2
	ONRS-0522	5 ml	22	1 inch	1:0.2
	ONRS-05221/2	5 ml	22	1½ inch	1:0.2
	ONRS-1021	10 ml	21	1 inch	1:0.4
	ONRS-10211/2	10 ml	21	1½ inch	1:0.4

- ONCE Concept syringe could change size syringe and needle that depend on your requirement.

ONCE™ Hypodermic needle

Standard range: 32G-16G, incorporate 304 s/s needle, PE needle cover with different color hub for easy recognition of size, Pyrogen Free Sterilized by EO Gas, Disposable use, Valid for 5 years.



All ONCE hypodermic needle are Triple Bevel Cut Sharpness it reduce size of puncture and pain when penetrated to the human skin.



ITEM NO.		COLOR		LENGTH							ITEM NO.
OD. (mm)	GAUGE	COLOR & Code		13mm	16mm	20mm	25mm	32mm	38mm	50mm	
				1/2"	5/8"	3/4"	1"	1 1/4"	1 1/2"	2"	
0.23	32G	Green	Green								ONHN-32G
0.26	31G	Orange	Orange								ONHN-31G
0.3	30G	Yellow	Yellow								ONHN-30G
0.33	29G	Red	Red								ONHN-29G
0.4	27G	Grey	Grey								ONHN-27G
0.45	26G	Brown	Brown								ONHN-26G
0.5	25G	Orange	Orange								ONHN-25G
0.55	24G	Purple	Purple								ONHN-24G
0.66	23G	Blue	Blue								ONHN-23G
0.7	22G	Black	Black								ONHN-22G
0.8	21G	Green	Green								ONHN-21G
0.9	20G	Yellow	Yellow								ONHN-20G
1.1	19G	Cream	Cream								ONHN-19G
1.2	18G	Pink	Pink								ONHN-18G
1.5	17G	White	White								ONHN-17G
1.65	16G	Lavender	Lavender								ONHN-16G

- Package: unit blister packing, bulk packing is optional
- Carton packing size: 543*380*528 mm, 8000PCS/CTN, 257CTN/20GP, 624CTN/40HQ

once
Intelligent Innovation

ONCE MEDICAL CO.,LTD
923/70P. SOI EAKACHAI30 EAKACHAI Rd. MAHACHAI
MUANG SAMUTSAKHON THAILAND 74000
TEL:+66(0) 34 422 449 FAX:+66(0) 34 413514 (AUTO)
WWW.ONCEMEDICALDEVICE.COM